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
FRAUNHOFER-INSTITUTE FOR APPLIED INFORMATION TECHNOLOGY FIT

Highlights

- Works in combination with fluorescence sensor technology
- Highly sensitive method
- Miniaturized, transportable system
- Reproducible handling of very small samples (<150 nl)
- Automated processes (enzyme assays, tissue staining using antibodies)

Fraunhofer Institute for Applied Information Technology FIT

Schloss Birlinghoven 53757 Sankt Augustin Germany

Biomolecular Optical Systems
Prof. Dr. Harald Mathis
Phone +49 2241 14-3514
harald.mathis@fit.fraunhofer.de

Contact:

Lorenz Sparrenberg Phone +49 2241 14-3752 lorenz.sparrenberg@fit.fraunhofer.de

www.fit.fraunhofer.de

POINT-OF-CARE DIAGNOSTICS

In drug research it is essential to understand the interactions of a pharmaceutical agent with endogenous substances. Until now this interaction could be observed only after the reaction – how the interaction evolved in detail remained hidden.

The electrowetting-on-dielectric system EWOD BioPro makes it possible to observe chemical reactions from start to finish. This allows measuring the kinetics of an interaction. In electrowetting, electric fields act on the electric charges in droplets, producing precise movements of individual droplets. The system works as follows:

- Individual droplets of a sample are created from a reservoir.
- A sample droplet is transported to a detection spot and merged with a solution droplet.
- The reaction is measured using appropriate sensors.
- The measurement data are evaluated and visualized with statistical methods.

Our EWOD BioPro system is particularly well suited for point-of-care diagnostics carried out not in a centralized lab but next to the patient's bed or in a doctor's office.

Implemented as a lab-on-a-chip system, it can analyze patient samples quickly and easily.

Application fields

- Point-of-care diagnostics
 - Examination of patient samples, e.g. sudor, saliva or urine, using enzymatic assays
 - Analysis of DNA samples by means of biological probes
- Parallelization and miniaturization of lab analyses
 - Reduced consumption of chemicals and sample materials
 - Faster examination results
- Basic research
 - Detailed understanding of biological processes